



Community Tissue Services™
Commitment • Technology • Safety

Commitment to serve our communities is of vital importance. Each day donors, their families, recipients and medical professionals are guided by our experienced tissue bank team members.

Technology drives advancement. As a result comes STERILE R dental allografts. Following a rigorous screening process, all dental allografts undergo the Allowash® process. This patented process facilitates the removal of cellular elements from musculoskeletal tissue while maintaining structural integrity. The Allowash® process utilizes both mechanical and chemical methodologies to reduce the potential of disease transmission.

All dental allografts are treated with irradiation in their final packaging.

Safety will not be compromised. All donors are recovered aseptically. Before dental allografts are ever released, all donors must pass a stringent screening process. This process meets or exceeds the Food and Drug Administration (FDA) regulations and the American Association of Tissue Banks (AATB) Standards.

衛署核准字號：
衛署醫字第0980214181號

Resorbable Membrane - Natural Barrier Grafts

STERILE R




Pericardium and Fascia Lata

Natural biocompatible resorbable membranes such as Pericardium and Fascia Lata are known for their ease of handling and protective capabilities. Membranes inhibit epithelial down growth and infiltration of connective tissue fibroblasts. Excellent for GTR & GBR procedures. Coforms to the graft site without memory.

Natural Barrier Grafts

Item No.	Description
DF152	1.5 X 2.0cm Fascia membrane
DF202	2.0 X 2.0cm Fascia membrane
DF203	2.0 X 3.0cm Fascia membrane
DF304	3.0 X 4.0cm Fascia membrane
DP152	1.5 X 2.0cm Pericardium membrane
DP202	2.0 X 2.0cm Pericardium membrane
DP203	2.0 X 3.0cm Pericardium membrane
DP304	3.0 X 4.0cm Pericardium membrane



Specialty Grafts



Mandible

Consider a natural solution in jaw replacement. Our experienced teams recover mandibles for this special need. Half Mandibles are available frozen and freeze dried. Whole Mandibles are only available frozen.



Costal Cartilage

Commonly used for facial and ridge defects, Costal Cartilage is available saline suspended and may be stored at room temperature.



RAY GARD

總代理：雷家有限公司

地址：407 台中市 中港路2段13-13號

電話：0800-324-925 • 02-8512-1511

傳真：02-8512-4300